

GREEN KERATIN EXTRACTION FROM DUCK FEATHERS USING A DEEP EUTECTIC SOLVENT

Sandra Alvarez, Eduardo Robles

▶ To cite this version:

Sandra Alvarez, Eduardo Robles. GREEN KERATIN EXTRACTION FROM DUCK FEATHERS USING A DEEP EUTECTIC SOLVENT. 3rd International Meeting on Deep Eutectic Systems, Jun 2023, Lisboa, Portugal. hal-04295340

HAL Id: hal-04295340 https://univ-pau.hal.science/hal-04295340

Submitted on 20 Nov 2023

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.
















GREEN KERATIN EXTRACTION FROM DUCK FEATHERS USING A DEEP EUTECTIC SOLVENT

Sandra Alvarez a*, Eduardo Robles a

^a University of Pau and the Adour Region, E2S UPPA, CNRS, Institute of Analytical and Physicochemical Sciences for the Environment and Material (IPREM-UMR 5254), 403 rue de Saint Pierre, 40000 Mont de Marsan, France


MATERIALS & METHODS

BACKGROUND



South the West of France, thousand of tons of duck feathers are produced by the food industry each year.

These feathers represent environmental and health risk because of their disposal in landfill sites or their incineration. However, a feather is composed of keratin: insoluble cysteinerich protein. It's high stability is ensured intramolecular / intermolecular disulfide bonds and hydrogen bonds. The main extraction methods employed are environmentally harmful and the demand to develop more eco-friendly processes is increasing. This work consists of determining an extraction process using a deep eutectic solvent and determining the influence of the extraction parameters on the yield and properties of the regenerated keratin.



DISCUSSION


- > A minimum of 130°C is required for extraction, increasing the temperature above 140°C reduces the overall yield. The shorter the dissolution time, the higher the yield. Keratin degradation is therefore minimized.
- > The optimum feather : solvent (F:S) ratio is 0.1. The limit results from the viscosity of the solution.
- Best thermal properties for lowest extraction temperature degradation from 240°C and highest remaining mass value. Thermal properties decrease overall with increasing temperature.
- Similar chemical composition for each temperature with characteristic protein bands evident.
- Use of non-toxic products for keratin environmentally-friendly extraction process has been optimized.




Keratin extraction and regeneration process Step 1 Step 2 Addition of Dissolution with Solvent elaboration Keratin regeneration grinded feathers Choline chloride – urea (1:2) different parameters



T= 60°C during 45 min



120°C < Temperature < 170 °C 90 min < time < 240 min



0,02 < F:S ratio < 0,15



Filtration

Step 3



Freeze-drying Step 5



Keratin neutralization by dialysis Step 4

RESULTS

Regenerated keratin





Figure 1. Yield as a function of temperature and extraction time, as well as the feather:solvent ratio (F:S).



Temperature (°C) Figure 3. Thermogravimetric analyses (TGA) of regenerated keratin depending on the extraction temperature.




Figure 2. X-ray diffraction (XRD) patterns of regenerated keratin.




Figure 4. Fourier-transform infrared (FTIR) spectra of regenerated keratin depending on the extraction temperature.

ACKNOWLEDGMENTS

The authors gratefully acknowledge the financial support of the Landes Departmental Council (CD40) and the tenure track position "Bois: Biobased materials", part of E2S UPPA supported by the "Investissements d'Avenir" French program managed by ANR (ANR-16-IDEX-0002).

